

Customer Repair Report

Model	Serial Number	Manufacture Dat	Manufacture Date		IRIDEX Service Contact	Completion Date	
OTHER	1333607			02715	Mark w		
Customer Name		Customer Number	Repair Type				
St. Marks Hospital		23375		☐ Warranty ☐ Upgrade ☐ Other			

Customer Description of Problem

BROKEN KNOBS

Technician Product Evaluate:

Cannot test due to Twisted and broken assembly fiber optic, wrong burnt halogen, broken control stem, extension cable read low power 1.20 KLx spec is 2.50 KLx white light shows dark shadow, loose viewing assembly from T-slot adapter, loose buld housing holder and will need alignment.

Service Performed on Unit: Waiting for customer approval

ITEMIZED REPAIR CHARGES

Part Number	Price/Unit	Quantity	Description	Total Price	
88080		1	Jumper Assy Mulitlam LIO+		984.00
10784		1	LIO Halogen Bulb-10W/6V Heine: X-004.88.068	\$	135.00
11481		1	DUST CAP W/LANYARD DART: IRIS DWG:	\$	3.00
87603		1	LIO EXTENSION CABLE(CS)	\$	812.00
TECHNICIAN HQ		5	Technician Labor	\$	1,425.00
BILLABLE FREIGHT		1	BILLABLE FREIGHT	\$	81.69

TOTAL CHARGES

3,440.69

Repaired products may contain reconditioned parts, but those are subject to the same specification and quality control standards applied to new materials and are warranted as such. The labor and parts are guaranteed for 90 days after the ship date of the product. Items not covered in the warranty are stated in the IRIDEX Standard Warranty Statement. Your signature below acknowledges the responsibility for payment of the listed charges. Freight and sales tax will be added where applicable.

Signature: Date: