

Customer Repair Report

Model	Serial Number		Manufacture Date		RMA Number	IRIDEX Service Contact	Completion Date	
OTHER	1265998				02721	Mark W.		
Customer Name			mer Number	Repair Type				
Intermountain Medical			21744	☐ Warranty ☐ Upgrade ☐ Other				

Customer Description of Problem

Preventative Maintenance

Technician Product Evaluate:

Cannote test due to Twisted and burnt assembly fiber optic, broken control stem, burnt halogen bulb, loose viewing assembly from T-slot adapter and will need alignment.

Service Performed on Unit: Waiting for customer approval

ITEMIZED REPAIR CHARGES

Part Number Price/Unit Q		Quantity	Description		Total Price	
10784		1	LIO Halogen Bulb-10W/6V Heine: X-004.88.068	\$	135.00	
11481		1	DUST CAP W/LANYARD DART: IRIS DWG:	\$	3.00	
12685		1	CONTROL STEM KNOB L & R VIEWING P/N: 11.65 38	\$	110.00	
12687		1	BEAM CONTROL STEMS HEINE SUPPORT KNOBS P/N: 11.65 15	\$	125.00	
65867		1	ASSY JUMPER LIO RFID	\$	984.00	
TECHNICIAN HQ		4.5	Technician Labor	\$	1,282.50	
BILLABLE FREIGHT	•	1	BILLABLE FREIGHT	\$	81.69	

TOTAL CHARGES

2,721.19

Repaired products may contain reconditioned parts, but those are subject to the same specification and quality control standards applied to new materials and are warranted as such. The labor and parts are guaranteed for 90 days after the ship date of the product. Items not covered in the warranty are stated in the IRIDEX Standard Warranty Statement. Your signature below acknowledges the responsibility for payment of the listed charges. Freight and sales tax will be added where applicable.

Signature: Date: