

Customer Repair Report

Model	Serial Number		Manufacture Date		RMA Number	IRIDEX Service Contact	Completion Date		
OTHER	L2788				02814	Mark W			
Customer Name			mer Number	Repair Type					
Children's Medical Center of Dallas			18136	Warranty Mon-warranty Upgrade Other					

Customer Description of Problem

unit overheating when using

Technician Product Evaluate:

Cannot test due to Twisted and burnt assembly fiber optic, wrong halogen bulb, bent control stem, loose viewing assembly from T-slot adapter and will need alignment. Not problem found over heater, LIO use halogen bulb is normal to be warm. Waiting for 88080 assembly fiber optic.

Service Performed on Unit: Waiting for customer approval

ITEMIZED REPAIR CHARGES

Part Number	· Price/Unit Quantity		Description		Total Price	
88080		1	Jumper Assy Mulitlam LIO+	\$	984.00	
10784		1	LIO Halogen Bulb-10W/6V Heine: X-004.88.068	\$	135.00	
11481		1	DUST CAP W/LANYARD DART: IRIS DWG:	\$	3.00	
12685		1	CONTROL STEM KNOB L & R VIEWING P/N: 11.65 38	\$	110.00	
12687		1	BEAM CONTROL STEMS HEINE SUPPORT KNOBS P/N: 11.65 15	\$	125.00	
TECHNICIAN HQ		4.5	Technician Labor	\$	1,282.50	
BILLABLE FREIGHT		1	BILLABLE FREIGHT	\$	191.60	

TOTAL CHARGES

2,831.10

Repaired products may contain reconditioned parts, but those are subject to the same specification and quality control standards applied to new materials and are warranted as such. The labor and parts are guaranteed for 90 days after the ship date of the product. Items not covered in the warranty are stated in the IRIDEX Standard Warranty Statement. Your signature below acknowledges the responsibility for payment of the listed charges. Freight and sales tax will be added where applicable.

Signature: Date: