

Customer Repair Report

Model	Serial Number		Manufacture Date		RMA Number	IRIDEX Service Contact	Completion Date	
OTHER	1362284				03050	Mark W		
Customer Name		Custon	omer Number Repair Type					
UT Medicine San Antonio			15115	☐ Warranty ☑ Upgrade ☐ Other				

Customer Description of Problem

UNABLE TO CENTER THE IMAGE

Technician Product Evaluate:

Cannot test due to viewing assembly dislocated from T-slot adapter, bent control control stem, burnt halogen bulb, T-slot adapter fractured, broken headband overband Omega 500 and will need alignment.

Service Performed on Unit: Waiting for customer approval

ITEMIZED REPAIR CHARGES

Part Number	Price/Unit	Quantity	Description		Total Price	
10784		1	LIO Halogen Bulb-10W/6V Heine: X-004.88.068	\$	135.00	
11481		1	DUST CAP W/LANYARD DART: IRIS DWG:	\$	3.00	
12685		1	CONTROL STEM KNOB L & R VIEWING P/N: 11.65 38	\$	110.00	
12687		1	BEAM CONTROL STEMS HEINE SUPPORT KNOBS P/N: 11.65 15	\$	125.00	
32185		1	OVERBAND OMEGA 500 MODIFIED MAKE FROM 32179	\$	177.96	
30800		1	T-SLOT ML LIO IRIDEX: 30800	\$	31.00	
TECHNICIAN HQ		4.5	Technician Labor	\$	1,282.50	
BILLABLE FREIGHT		1	BILLABLE FREIGHT	\$	191.60	

TOTAL CHARGES \$

2,056.06

Repaired products may contain reconditioned parts, but those are subject to the same specification and quality control standards applied to new materials and are warranted as such. The labor and parts are guaranteed for 90 days after the ship date of the product. Items not covered in the warranty are stated in the IRIDEX Standard Warranty Statement. Your signature below acknowledges the responsibility for payment of the listed charges. Freight and sales tax will be added where applicable.

Print Name:	
Signature:	Date: