

Customer Repair Report

Model	Serial Number		Manufacture Date		RMA Number	IRIDEX Service Contact	Completion Date	
OTHER	1171644				03149	Mark W		
Customer Name			mer Number	Repair Type				
Alta Bates Medical Center			12833	Warranty Mon-warranty Upgrade Other				

Customer Description of Problem

The mirror will not hold is position

Technician Product Evaluate:

Cannot test due to Twisted and damaged extension cable, wrong and buert halogen bulb, broken control stem, loose viewing assembly from T-slot adapter, scratched assembly fiber optic and will need alignment.

Service Performed on Unit: Waiting for customer approval

ITEMIZED REPAIR CHARGES

Part Number	Part Number Price/Unit Quantit		Description	Total Price	
10784		1	LIO Halogen Bulb-10W/6V Heine: X-004.88.068	\$	135.00
11481		1	DUST CAP W/LANYARD DART: IRIS DWG:	\$	3.00
87603		1	LIO EXTENSION CABLE(CS)	\$	812.00
12687		1	BEAM CONTROL STEMS HEINE SUPPORT KNOBS P/N: 11.65 15	\$	125.00
12685		1	CONTROL STEM KNOB L & R VIEWING P/N: 11.65 38	\$	110.00
TECHNICIAN HQ		5	Technician Labor	\$	1,425.00
BILLABLE FREIGHT		1	BILLABLE FREIGHT	\$	53.60

TOTAL CHARGES

2,663.60

Repaired products may contain reconditioned parts, but those are subject to the same specification and quality control standards applied to new materials and are warranted as such. The labor and parts are guaranteed for 90 days after the ship date of the product. Items not covered in the warranty are stated in the IRIDEX Standard Warranty Statement. Your signature below acknowledges the responsibility for payment of the listed charges. Freight and sales tax will be added where applicable.

Signature: Date: